

iCUP Rx

iCup Rx



DRUG TESTING KITS

- Screen for 5 commonly abused prescription drugs
- Detects Opiates buprenorphine at low cut-off levels
- CLIA-Waived
- Results in only 5 minutes
- Peel label to reveal results
- Test card built into collection cup
- Simple to use, no need to tilt cup
- Specimen validity test and temperature strip included

ALCOPRO
DRUG & ALCOHOL TESTING PRODUCTS

1-800-227-9890
www.alcopro.com

iCup Rx



iCup Rx 5 Panel Drug Test

iCup Rx 5 panel drug test is an urine drug test that targets 5 of the most commonly abused prescription drugs and detects Opiates at low-levels. Great for clinicians, drug addiction professionals and probationary programs, the iCup Rx is the most reliable screening test where compliance must be frequently monitored. Simple to use, there's no need to tilt the cup to activate the test and results are ready in 5 minutes.

Easy to Use

The iCup Rx is CLIA waived and uses an integrated cup design that builds a drug detection card securely into a specimen container to minimize contact with the specimen. The donor provides a urine specimen directly into the iCup and attaches the leak-proof, screw-on lid. Specimen Validity Test and a built in temperature strip helps monitor specimen tampering.

Easy Laboratory Confirmation

Send the leak proof iCup Rx container with the urine specimen to a laboratory for confirmation testing of non-negative results. This avoids the need to split or pour the specimen into another container, saving time and avoiding additional handling of the specimen. Using the security seals included with the iCup creates a strong chain of custody that documents from the point of collection to sending the specimen to the laboratory for confirmation testing.

EASE OF USE



Step 1

Temperature strip assures valid specimen.



Step 2

Peel back label to view vivid, dark result lines.



Step 3

Leak-proof lid allows operator to lay specimen on it's side for photocopying results.

iCup Rx

Item	Description	1 - 24	25+
8iRx	Bup/Bzo/Mtd/Opi/Oxy	\$ 6.50 ea	\$5.75 ea

Specifications	
Test Medium	Urine
Format	Integrated Split specimen cup with temperature strip
Cut-off Levels	
Bup (buprenorphine)	10 ng/ml
Bzo (benzodiazepine)	300 ng/ml
MTD (methadone)	300 ng/ml
OPI (opiate)	300 ng/ml
Oxy (oxycodone)	100 ng/ml

For Information & Orders Call: **1-800-227-9890**
 Order 24 Hours A Day At: **www.alcopro.com**
 Fax Orders To: **1-800-655-5834**

ALCOPRO
 DRUG & ALCOHOL TESTING PRODUCTS